

## [Circle Hereditary Cancer Genetic Test White Paper – last updated 1.11.2018](#)

### Executive Summary

Prenetics has developed a targeted next generation sequencing (NGS) based test for hereditary cancer. We developed this test using Illumina's NGS platform by analyzing 35 hereditary cancer genes. This test is designed to identify cancer-causing mutations for 8 different cancers - breast, ovarian, prostate, uterine, colorectal, pancreatic, stomach cancers and melanoma. In this test, we sequence the protein-coding regions of the 35 genes. We have also developed a bioinformatics pipeline that utilizes standard practices in the field. The test has a high analytical validity for the detection of single nucleotide variants and small insertions and deletions (indels). Upon testing 2976 variants across 58 samples, the analytical accuracy of the test is 99.8%.

### Introduction

Understanding the genetic basis of cancer risk is a major international endeavor that started in 1980s. The identification of mutations in the genes BRCA1 and BRCA2 in the 1990s was a major turning point in the acceleration of research in hereditary cancer. Advent of next generation sequencing (NGS) in late 2000's has further accelerated the discovery of many cancer susceptibility genes. More than 100 cancer susceptibility genes that are associated with cancer risk have been discovered (Stadler et al 2014) and more genes are being discovered continuously (Sokolenko et al 2015).

While discovery of such new cancer genes is a critical first step, additional clinical studies are needed to independently validate these discoveries. However, development of such clinical evidence for the new genes typically lags behind the discovery of these genes. Based on our assessment of strength of available clinical evidence, we have generated a hereditary cancer panel including 35 genes.

In order to assess strength of clinical evidence for the genes included in our cancer panel, we assess the following sources: a) Testing of the cancer genes as recommended by professional guidelines (e.g., from National Comprehensive Cancer Network (NCCN), American Society of Clinical Oncology (ASCO)), b) Genes found to be valid cancer genes by systematic reviews, c) Genes in which pathogenic mutations have been reported by multiple research studies or reputed resources such as ClinGen.

Based on the above criteria, we selected 35 genes related to 8 cancers - breast, ovarian, prostate, uterine, colorectal, pancreatic, stomach cancers and melanoma. In order to analytically validate the panel, we tested our panel using reference samples. Below, we provide the detailed methods and results of our validation process.

### Materials & Methods

Prenetics has developed a next generation sequencing (NGS) based hereditary cancer test and has followed a process of established NGS laboratory protocols and custom bioinformatics analysis to generate the test results. Below, we briefly describe our procedures: Genomic DNA is extracted from saliva sample by standard extraction method. Extracted DNA is fragmented and prepared into NGS libraries compatible with the Illumina platform. The Nimblegen SeqCap Target Enrichment technology is applied to enrich 35 genes selected based on their association with different cancers. The target enriched DNA fragments are analyzed by utilizing 2X150 paired-end

sequencing method with an Illumina MiSeq Sequencer with at least 50X average coverage to ensure data quality. Quality control procedures are in-place along the process to ensure sample identification, high quality of DNA isolation, library preparation, target capture and sequencing. Each sequencing test includes one fully characterized positive control. The standard bioinformatics analysis pipeline is used to analyze sequencing data with reference to genome GRCh37 (hg19) for sequence alignment. Our pipeline includes BWA alignment (0.7.15) and GATK (4.0) best practices. Low quality and duplicated reads are removed and variants are detected by following a well-established bioinformatics variants calling pipeline. Variants are classified by following the standards and guidelines for sequence variants interpretation of the American College of Medical genetics and Genomics (ACMG). Variant annotation is performed by using the ANNOVAR framework. Our pipeline is built using a Python library called Snakemake (Koster et al 2012). Snakemake is a workflow management system that is designed to create reproducible and scalable data analyses.

Variants are classified into five categories, including the pathogenic, likely pathogenic, variant of uncertain significance (VUS), likely benign, and benign. Pathogenic and likely pathogenic variants identified by NGS undergo additional confirmatory test using Sanger sequencing before getting reported. All VUS results are reviewed bi-annually for updates in the scientific literatures and notifications will be sent if there are changes to the classification of reported variants.

Our hereditary cancer test is designed to assess clinically relevant mutations in 35 genes associated with hereditary risk for breast, ovarian, colorectal, melanoma, pancreatic, prostate, uterine and stomach cancer (detailed in Supplement Table 1) through the detection of single nucleotide variants and small insertions and deletions (Indel; up to 15-25bp) located in the DNA coding sequences, nearby flanking regions (+/- 20 bp flanking of each exon) and known splice regions in the genes targeted by the cancer panel. These genes include *APC*, *ATM*, *BAP1*, *BARD1*, *BMPR1A*, *BRCA1*, *BRCA2*, *BRIP1*, *CDH1*, *CDK4*, *CDKN2A*, *CHEK2*, *EPCAM*, *GREM1*, *MEN1*, *MITF*, *MLH1*, *MRE11A*, *MSH2*, *MSH6*, *MUTYH*, *NBN*, *PALB2*, *PMS1*, *PMS2*, *POLD1*, *POLE*, *PTEN*, *RAD50*, *RAD51C*, *RAD51D*, *SMAD4*, *STK11*, *TP53*, and *XRCC2*. These genes are assessed for variants within all coding exons. This panel is not intended to analyze the following types of variations: gross rearrangements and deep intronic variation, Alu element insertions, certain homopolymers such as those in *PMS1* and other unknown abnormalities. The pattern of variant types varies with the genes and this test detects a high but variable percentage of known and unknown variants of the classes stated. Novel sequence changes in the promoter regions and other non-coding regions will not be detected by this test.

Our validation strategy is governed by the guidelines for NGS from the ACMG (Rehm et al 2013) and the College of American Pathologist (CAP). The validation studies are designed to provide a good representation of possible variant types across the 35 gene in our cancer panel. These groups included well-characterized cell lines and/or DNA specimens, previously extracted from blood from individuals whose genetic variants had been previously characterized by the 1000 Genome Project (The 1000 Genome Project Consortium; Sudmant et al 2015)

## Validation Results

### 1) Descriptive statistics of variants from reference materials

A total of 58 reference materials are selected based on the availability of high quality variant calls and recommendations by the National Institute of Standards and Technology (NIST) (Zook et al 2014). Variant calls in these reference materials were compared to 1000 Genome Project and

Illumina Platinum Genomes (The 1000 Genome Project Consortium; Sudmant et al 2015; Eberie et al 2017).

Table 1. Descriptive variant statistics from reference materials, organized by individual samples.

Specimen	Number of variants		Total	Specimen	Number of variants		Total	Specimen	Number of variants		Total
	SNVs	Indels			SNVs	Indels			SNVs	Indels	
HG00625	48	1	49	HG00701	51	1	52	NA18561	54	2	56
HG00626	54	0	54	HG00704	59	1	60	NA18562	57	0	57
HG00628	46	1	47	HG00705	49	1	50	NA18563	40	0	40
HG00629	48	0	48	HG00707	51	1	52	NA18572	57	0	57
HG00634	61	1	62	HG00708	61	0	61	NA18603	53	2	55
HG00650	49	0	49	HG02521	56	0	56	NA18605	57	0	57
HG00651	58	2	60	HG02522	59	1	60	NA18606	46	1	47
HG00653	35	1	36	NA12878	53	1	54	NA18608	45	1	46
HG00654	57	0	57	NA12889	53	1	54	NA18609	56	1	57
HG00656	58	1	59	NA12890	39	1	40	NA18611	55	0	55
HG00657	39	0	39	NA18530	50	0	50	NA18612	46	0	46
HG00662	49	1	50	NA18534	54	1	55	NA18613	53	1	54
HG00663	55	1	56	NA18536	36	1	37	NA18620	52	0	52
HG00671	55	1	56	NA18543	40	1	41	NA18621	46	1	47
HG00672	36	1	37	NA18544	48	1	49	NA18622	43	1	44
HG00683	57	0	57	NA18546	46	0	46	NA18623	58	0	58
HG00684	65	1	66	NA18548	60	0	60	NA18624	50	1	51
HG00689	46	1	47	NA18549	50	1	51	NA18629	50	0	50
HG00690	48	1	49	NA18557	54	1	55	Subtotal	918	11	929
HG00698	43	1	44	NA18558	41	1	42	<b>TOTAL</b>	<b>2935</b>	<b>41</b>	<b>2976</b>
Subtotal	1007	15	1022	Subtotal	1010	15	1025				

## 2) Independent confirmation of variants

The 35-gene cancer genetic test was validated using 58 specimens selected from the 1000 genome project (Coriell Institute for Medical Research). All functional types of variants are included for this analysis. The results of our cancer panel was compared against the previously identified variants in 1000 genomes project which allowed the assessment of test accuracy. Low-confidence variants in 1000 genomes project were independently confirmed by Sanger sequencing. Table 2 below provides the results of this analysis organized by genes.

Table 2. Overview of results for 2976 variants.

Gene	Total variants	True Positives	False Positives	False Negatives
<i>APC</i>	395	395	0	0
<i>ATM</i>	73	73	0	0
<i>BAP1</i>	11	11	0	0
<i>BARD1</i>	178	178	0	0
<i>BMPR1A</i>	52	52	0	0
<i>BRCA1</i>	229	229	0	0
<i>BRCA2</i>	359	359	0	1
<i>BRIP1</i>	174	174	0	0
<i>CDH1</i>	69	69	0	0
<i>CDKN2A</i>	76	76	0	0
<i>CHEK2</i>	12	12	0	0
<i>EPCAM</i>	73	73	0	0
<i>GREM1</i>	24	24	0	0
<i>MEN1</i>	154	154	0	0
<i>MITF</i>	1	1	0	0
<i>MLH1</i>	17	17	0	0
<i>MRE11A</i>	90	90	0	0
<i>MSH2</i>	61	61	0	0
<i>MSH6</i>	54	54	0	0
<i>MUTYH</i>	43	43	0	0
<i>NBN</i>	180	180	0	0
<i>PALB2</i>	26	26	0	0
<i>PMS1</i>	31	31	0	0
<i>PMS2</i>	232	232	0	4*
<i>POLD1</i>	51	51	0	0
<i>POLE</i>	211	211	0	0
<i>PTEN</i>	6	6	0	0
<i>RAD50</i>	2	2	0	0

<i>RAD51C</i>	7	7	0	0
<i>RAD51D</i>	29	29	0	0
<i>STK11</i>	7	7	0	0
<i>TP53</i>	48	48	0	0
<i>XRCC2</i>	1	1	0	0
<b>TOTAL</b>	<b>2976</b>	<b>2976</b>	<b>0</b>	<b>5</b>

\*4 variants at the position Chr7:6013153 in PMS2 is a false negative. But, this is in a pseudogene region and is also a benign variant.

## Interpretation and Conclusion

Our cancer panel achieves an analytical accuracy of 99.8% for single nucleotide variants and small indels (<25bp) based on the above validation study results. Sensitivity to detect insertions and deletions larger than 25bp but smaller than a full exon may be marginally reduced. Overall, the test has a high analytical validity for the detection of single nucleotide variants and small insertions and deletions.

## References

- 1000 Genomes Project Consortium. A global reference for human genetic variation. *Nature*. 2015 Oct 1;526(7571):68-74.
- Aloraifi F, McDevitt T, Martiniano R, McGreevy J, McLaughlin R, Egan CM, Cody N, Meany M, Kenny E, Green AJ, Bradley DG. Detection of novel germline mutations for breast cancer in non-BRCA1/2 families. *The FEBS journal*. 2015 Sep 1;282(17):3424-37.
- Al-Tassan N, Chmiel NH, Maynard J, Fleming N, Livingston AL, Williams GT, Hodges AK, Davies DR, David SS, Sampson JR, Cheadle JP. Inherited variants of MYH associated with somatic G: C→ T: A mutations in colorectal tumors. *Nature genetics*. 2002 Feb 1;30(2):227-32.
- BRCA SG. A strong candidate for the breast and ovarian cancer susceptibility gene BRCA1. *Science*. 1994 Oct;266:7.
- Bronner CE, Baker SM, Morrison PT, Warren G, Smith LG, Lescoe MK, Kane M, Earabino C, Lipford J, Lindblom A, Tannergrnd P. Mutation in the DNA mismatch repair gene homologue hMLH 1 is associated with hereditary non-polyposis colon cancer. *Nature*. 1994 Mar 17;368(6468):258-61.
- Castéra L, Krieger S, Rousselin A, Legros A, Baumann JJ, Bruet O, Brault B, Fouillet R, Goardon N, Letac O, Baert-Desurmont S. Next-generation sequencing for the diagnosis of hereditary breast and ovarian cancer using genomic capture targeting multiple candidate genes. *European Journal of Human Genetics*. 2014 Nov 1;22(11):1305-13.
- Chandrasekharappa SC, Guru SC, Manickam P, Olufemi SE, Collins FS, Emmert-Buck MR, Debelenko

- LV, Zhuang Z, Lubensky IA, Liotta LA, Crabtree JS. Positional cloning of the gene for multiple endocrine neoplasia-type 1. *Science*. 1997 Apr 18;276(5311):404-7.
8. Eberle MA, Fritzilas E, Krusche P, Källberg M, Moore BL, Bekritsky MA, Iqbal Z, Chuang HY, Humphray SJ, Halpern AL, Kruglyak S. A reference data set of 5.4 million phased human variants validated by genetic inheritance from sequencing a three-generation 17-member pedigree. *Genome research*. 2017 Jan 1;27(1):157-64.
  9. Erkkö H, Xia B, Nikkilä J, Schleutker J, Syrjäkoski K, Mannermaa A, Kallioniemi A, Pylkäs K, Karppinen SM, Rapakko K, Miron A. A recurrent mutation in PALB2 in Finnish cancer families. *Nature*. 2007 Mar 15;446(7133):316-9.
  10. Fishel R, Lescoe MK, Rao MR, Copeland NG, Jenkins NA, Garber J, Kane M, Kolodner R. The human mutator gene homolog MSH2 and its association with hereditary nonpolyposis colon cancer. *Cell*. 1993 Dec 3;75(5):1027-38.
  11. Groden J, Thliveris A, Samowitz W, Carlson M, Gelbert L, Albertsen H, Joslyn G, Stevens J, Spirio L, Robertson M, Sargeant L. Identification and characterization of the familial adenomatous polyposis coli gene. *Cell*. 1991 Aug 9;66(3):589-600.
  12. Guilford P, Hopkins J, Harraway J, McLeod M, McLeod N, Harawira P, Taite H, Scoular R, Miller A, Reeve AE. E-cadherin germline mutations in familial gastric cancer. *Nature*. 1998 Mar 26;392(6674):402-5.
  13. Hemminki A, Markie D, Tomlinson I, Avizienyte E, Roth S, Loukola A, Bignell G, Warren W, Aminoff M, Höglund P, Järvinen H. A serine/threonine kinase gene defective in Peutz–Jeghers syndrome. *Nature*. 1998 Jan 8;391(6663):184-7.
  14. Howe JR, Bair JL, Sayed MG, Anderson ME, Mitros FA, Petersen GM, Velculescu VE, Traverso G, Vogelstein B. Germline mutations of the gene encoding bone morphogenetic protein receptor 1A in juvenile polyposis. *Nature genetics*. 2001 Jun 1;28(2):184-7.
  15. Howe JR, Roth S, Ringold JC, Summers RW, Järvinen HJ, Sistonen P, Tomlinson IP, Houlston RS, Bevan S, Mitros FA, Stone EM. Mutations in the SMAD4/DPC4 gene in juvenile polyposis. *Science*. 1998 May 15;280(5366):1086-8.
  16. Hussussian CJ, Struewing JP, Goldstein AM, Higgins PA, Ally DS, Sheahan MD, Clark WH, Tucker MA, Dracopoli NC. Germline p16 mutations in familial melanoma. *Nature genetics*. 1994 Sep 1;8(1):15-21.
  17. Jenne DE, Reomann H, Nezu JI, Friedel W, Loff S, Jeschke R, Müller O, Back W, Zimmer M. Peutz–Jeghers syndrome is caused by mutations in a novel serine threonine kinase. *Nature genetics*. 1998 Jan 1;18(1):38-43.
  18. Joslyn G, Carlson M, Thliveris A, Albertsen H, Gelbert L, Samowitz W, Groden J, Stevens J, Spirio L, Robertson M, Sargeant L. Identification of deletion mutations and three new genes at the familial



polyposis locus. *Cell*. 1991 Aug 9;66(3):601-13.

19. Kamb A, Shattuck-Eidens D, Eeles R, Liu Q, Gruis NA, Ding W, Hussey C, Tran T, Miki Y, Weaver-Feldhaus J, McClure M. Analysis of the p16 gene (CDKN2) as a candidate for the chromosome 9p melanoma susceptibility locus. *Nature genetics*. 1994 Sep 1;8(1):22-6.
20. Köster, Johannes and Rahmann, Sven. "Snakemake - A scalable bioinformatics workflow engine". *Bioinformatics* 2012.
21. Leach FS, Nicolaides NC, Papadopoulos N, Liu B, Jen J, Parsons R, Peltomäki P, Sistonen P, Aaltonen LA, Nyström-Lahti M, Guan XY. Mutations of a mutS homolog in hereditary nonpolyposis colorectal cancer. *Cell*. 1993 Dec 17;75(6):1215-25.
22. Levitus M, Waisfisz Q, Godthelp BC, de Vries Y, Hussain S, Wiegant WW, Elghalbzouri-Maghrani E, Steltenpool J, Roomans MA, Pals G, Arwert F. The DNA helicase BRIP1 is defective in Fanconi anemia complementation group J. *Nature genetics*. 2005 Sep 1;37(9):934-5.
23. Levrán O, Attwooll C, Henry RT, Milton KL, Neveling K, Rio P, Batish SD, Kalb R, Velleuer E, Barral S, Ott J. The BRCA1-interacting helicase BRIP1 is deficient in Fanconi anemia. *Nature genetics*. 2005 Sep 1;37(9):931-3.
24. Liaw D, Marsh DJ, Li J, Dahia PL, Wang SI, Zheng Z, Bose S, Call KM, Tsou HC, Peacocke M, Eng C. Germline mutations of the PTEN gene in Cowden disease, an inherited breast and thyroid cancer syndrome. *Nature genetics*. 1997 May 1;16(1):64-7.
25. Litman R, Peng M, Jin Z, Zhang F, Zhang J, Powell S, Andreassen PR, Cantor SB. BACH1 is critical for homologous recombination and appears to be the Fanconi anemia gene product FANCI. *Cancer cell*. 2005 Sep 30;8(3):255-65.
26. Loveday C, Turnbull C, Ramsay E, Hughes D, Ruark E, Frankum JR, Bowden G, Kalmyrzaev B, Warren-Perry M, Snape K, Adlard JW. Germline mutations in RAD51D confer susceptibility to ovarian cancer. *Nature genetics*. 2011 Sep 1;43(9):879-82.
27. Malkin D, Li FP, Strong LC, Fraumeni Jr JF, Nelson CE, Kim DH, Kassel J, Gryka MA, Bischoff FZ, Tainsky MA, Friend SH. Germ line p53 mutations in a familial syndrome of breast cancer, sarcomas, and other neoplasms. *science*. 1990 Nov 30:1233-8.
28. Matsuura S, Tauchi H, Nakamura A, Kondo N, Sakamoto S, Endo S, Smeets D, Solder B, Belohradsky BH, Der Kaloustian VM, Oshimura M. Positional cloning of the gene for Nijmegen breakage syndrome. *Nature genetics*. 1998 Jun 1;19(2):179-81.
29. Meijers-Heijboer H, Van den Ouweland A, Klijn J, Wasielewski M, de Snoo A, Oldenburg R, Hollestelle A, Houben M, Crepin E, van Veghel-Plandsoen M, Elstrodt F. Low-penetrance susceptibility to breast cancer due to CHEK2\* 1100delC in noncarriers of BRCA1 or BRCA2 mutations. *Nature genetics*. 2002 May 1;31(1):55-9.

30. Meindl A, Hellebrand H, Wiek C, Erven V, Wappenschmidt B, Niederacher D, Freund M, Lichtner P, Hartmann L, Schaal H, Ramser J. Germline mutations in breast and ovarian cancer pedigrees establish RAD51C as a human cancer susceptibility gene. *Nature genetics*. 2010 May 1;42(5):410-4.
31. Miyaki M, Konishi M, Tanaka K, Kikuchi-Yanoshita R, Muraoka M, Yasuno M, Igari T, Koike M, Chiba M, Mori T. Germline mutation of MSH6 as the cause of hereditary nonpolyposis colorectal cancer. *Nature genetics*. 1997 Nov 1;17(3):271-2.
32. Nishisho I, Nakamura Y. Mutations of chromosome 5q21 genes in FAP and colorectal cancer patients. *Science*. 1991 Aug 9;253(5020):665.
33. Papadopoulos N, Nicolaidis NC, Wei YF, Ruben SM, Carter KC, Rosen CA, Haseltine WA, Fleischmann RD, Fraser CM, Adams MD, Venter JC. Mutation of a'mutL'homolog in hereditary colon cancer. *Science*. 1994 Mar 18;263(5153):1625-30.
34. Preisinger AC, Hedge PH, McKechnie DO, Finniear R, Markham A, Groffen J, Boguski M, Altschui S, Horii A, Ando H, Miyoshi Y. Identification of FAP locus genes from chromosome 5q21. *Science*. 1991 Aug 9;253(5020):661-5.
35. Rahman N. Realizing the promise of cancer predisposition genes. *Nature*. 2014 Jan 16;505(7483):302-8.
36. Rahman N, Seal S, Thompson D, Kelly P, Renwick A, Elliott A, Reid S, Spanova K, Barfoot R, Chagtai T, Jayatilake H. PALB2, which encodes a BRCA2-interacting protein, is a breast cancer susceptibility gene. *Nature genetics*. 2007 Feb 1;39(2):165-7.
37. Rehm HL, Bale SJ, Bayrak-Toydemir P, Berg JS, Brown KK, Deignan JL, Friez MJ, Funke BH, Hegde MR, Lyon E, Working Group of the American College of Medical Genetics. ACMG clinical laboratory standards for next-generation sequencing. *Genetics in medicine*. 2013 Jul 25;15(9):733-47.
38. Reid S, Schindler D, Hanenberg H, Barker K, Hanks S, Kalb R, Neveling K, Kelly P, Seal S, Freund M, Wurm M. Biallelic mutations in PALB2 cause Fanconi anemia subtype FA-N and predispose to childhood cancer. *Nature genetics*. 2007 Feb 1;39(2):162-4.
39. Ricker C, Culver JO, Lowstuter K, Sturgeon D, Sturgeon JD, Chanock CR, Gauderman WJ, McDonnell KJ, Idos GE, Gruber SB. Increased yield of actionable mutations using multi-gene panels to assess hereditary cancer susceptibility in an ethnically diverse clinical cohort. *Cancer genetics*. 2016 Apr 30;209(4):130-7.
40. Savitsky K, Bar-Shira A, Gilad S, Rotman G, Ziv Y, Vanagaite L, Tagle DA, Smith S, Uziel T, Sfez S, Ashkenazi M. A single ataxia telangiectasia gene with a product similar to PI-3 kinase. *Science*. 1995 Jun 23;1749-53.
41. Sokolenko AP, Suspitsin EN, Kuligina ES, Bizin IV, Frishman D, Imyanitov EN. Identification of novel hereditary cancer genes by whole exome sequencing. *Cancer letters*. 2015 Dec 28;369(2):274-88.



42. Stadler ZK, Schrader KA, Vijai J, Robson ME, Offit K. Cancer genomics and inherited risk. *Journal of Clinical Oncology*. 2014 Jan 21;32(7):687-98.
43. Sudmant PH, Rausch T, Gardner EJ, Handsaker RE, Abyzov A, Huddleston J, Zhang Y, Ye K, Jun G, Fritz MH, Konkel MK. An integrated map of structural variation in 2,504 human genomes. *Nature*. 2015 Oct 1;526(7571):75-81.
44. Susswein LR, Marshall ML, Nusbaum R, Postula KJ, Weissman SM, Yackowski L, Vaccari EM, Bissonnette J, Booker JK, Cremona ML, Gibellini F. Pathogenic and likely pathogenic variant prevalence among the first 10,000 patients referred for next-generation cancer panel testing. *Genetics in Medicine*. 2015 Dec 17;18(8):823-32.
45. Vahteristo P, Bartkova J, Eerola H, Syrjäkoski K, Ojala S, Kilpivaara O, Tamminen A, Kononen J, Aittomäki K, Heikkilä P, Holli K. A CHEK2 genetic variant contributing to a substantial fraction of familial breast cancer. *The American Journal of Human Genetics*. 2002 Aug 31;71(2):432-8.
46. Varon R, Vissinga C, Platzer M, Cerosaletti KM, Chrzanowska KH, Saar K, Beckmann G, Seemanová E, Cooper PR, Nowak NJ, Stumm M. Nibrin, a novel DNA double-strand break repair protein, is mutated in Nijmegen breakage syndrome. *Cell*. 1998 May 1;93(3):467-76.
47. Wiesner T, Obenauf AC, Murali R, Fried I, Griewank KG, Ulz P, Windpassinger C, Wackernagel W, Loy S, Wolf I, Viale A. Germline mutations in BAP1 predispose to melanocytic tumors. *Nature genetics*. 2011 Oct 1;43(10):1018-21.
48. Wooster R, Blignell G, Lancaster J, Surlfti S, Seal S. Identification of the breast cancer susceptibility gene BRCA2. *Nature*. 1995 Dec 28;378:21.
49. Xia B, Dorsman JC, Ameziane N, de Vries Y, Rooimans MA, Sheng Q, Pals G, Errami A, Gluckman E, Llera J, Wang W. Fanconi anemia is associated with a defect in the BRCA2 partner PALB2. *Nature genetics*. 2007 Feb 1;39(2):159-61.
50. Zook JM, Chapman B, Wang J, Mittelman D, Hofmann O, Hide W, Salit M. Integrating human sequence data sets provides a resource of benchmark SNP and indel genotype calls. *Nature biotechnology*. 2014 Mar 1;32(3):246-51.

## Supplement

Supplemental Table 1: Known associations between genes in our cancer panel and cancer type

Gene	Breast	Ovarian	Uterine	Colorectal	Melanoma	Pancreatic	Stomach	Prostate*
APC				•		•	•	
ATM	•			•		•	•	•
BAP1					•			
BARD1	•	•						
BMPR1A				•		•	•	
BRCA1	•	•				•		•
BRCA2	•	•			•	•		•
BRIP1	•	•						
CDH1	•			•			•	
CDK4					•			
CDKN2A					•	•		
CHEK2	•	•		•				•
EPCAM		•	•	•		•	•	
GREM1				•				
MEN1						•	•	
MITF					•			
MLH1		•	•	•		•	•	
MRE11A	•	•		•				
MSH2		•	•	•		•	•	
MSH6		•	•	•		•	•	
MUTYH	•		•	•			•	
NBN	•	•						•
PALB2	•	•				•		•
PMS1				•				
PMS2		•	•	•		•	•	•
POLD1				•				
POLE				•				
PTEN	•		•	•	•			
RAD50	•	•						
RAD51C	•	•						
RAD51D	•	•						
SMAD4				•		•	•	
STK11	•	•	•	•		•	•	
TP53	•	•	•	•	•	•	•	•
XRCC2	•	•		•				

This table depicts only primary cancer risks and does not specify other important gene-specific associated cancers.

\* Please note that research and screening guidelines for genes associated with hereditary prostate cancer are still in their early stages.